



Product Catalog 2014

New and Precise
High Sensitivity and Specificity
Stable Quality





Our Company

Healgen Scientific LLC is a professional bio-tech products manufacturer and distributor located in Houston, Texas. It has developed a range of diagnostic products, laboratory consumables as well as equipments. Our aim is to provide diagnostic laboratories and any other bio-tech R&D organizations with high quality products and services. It is always our premium to enable and inspire scientists and doctors to reach their R&D and clinical goals. Our highly trained scientists and engineers will fully understand your needs and requirements. By using our intellectual property rights, including: A Labelling Method of Pathogenic Cells Exogenetic Antigen (200610155201.3), Preparation of Indole and its derivatives (200710069695.9), A Target Gene Expression method with enhanced cell specificity (Patent No.: ZL1118525.2 International Patent classification No.: 12N 15/67), Healgen will always provide its core competency and help you develop your research ideas.

Our Products include:


Rapid Test, clinical Chemistry Reagent, Equipment, Second Antibodies, Pipette, FISH Kit, Laboratory Consumables, and some of them have got CE and FDA certificate.

Healgen Scientific has established one manufacturing facilities in China, named Zhejiang Orient Gene Biotech Co., Ltd. Zhejiang Orient Gene Biotech Co., Ltd is the main manufacturing facility of Healgen Scientific LLC. It covers 144 acres and plans to build 92,000 M2 Of GMP workshop, science lab and other workshops, accessory house, production and office house. Currently we have built 15,000 M2 GMP workshop and has over 300 employees.

Adopting the most advanced enzyme, cell and gene engineering technologies in the world, our diagnostic products not only have special advantages in diagnosing infectious diseases, but also have a good prospect in tumor diagnosing.

Moreover, most products of Zhejiang Orient Gene Biotech Co., Ltd are CE and SFDA certified. Products have been sold to USA, Europe, Mid-east, South Asia, Africa and South America.

We are also willing to take OEM orders and work as agent for other biotech companies. Supported by our professional personnel and good manufacturing practice, we will always offer the best quality and services to our customers.



3-5

Infectious Disease



6-10

Drugs of Abuse



13

Cardiac Marker



14

Urinalysis Strip

Rapid Test

Infectious Disease	3-5
Drugs of Abuse	6-10
Pregnancy & Fertility	11
Tumor Marker	12
Cardiac Marker	13
Urinalysis Strip	14
Veterinary Test	15

Infectious Disease



Product Name	Specimen	Catalog No.	Format	Sensitivity	Package	Qualification
Chlamydia Rapid Test	Swab/Urine	GCCHL-502a	Device	1.0E*7	20T/Box	CE
		GCCHL-504a	Card	1.0E*7		
Dengue Fever IgG/IgM antibody Rapid Test	S/P	GCDEN(ab)-301a	Strip	/	50T/Box	CE
		GCDEN(ab)-302a	Device	/	25T/Box	CE
		GCDEN(ab)-304a	Card	/		
	WB/S/P	GCDEN(ab)-401a	Strip	/	50T/Box	CE
		GCDEN(ab)-402a	Device	/	25T/Box	CE
		GCDEN(ab)-402b	Device	/	25T/Box	CE
		GCDEN(ab)-404a	Card	/		
		GCDEN(ab)-404b	Device	/	25T/Box	CE
Dengue Fever NS 1 antigen Rapid Test	S/P	GCDEN(NS)-302a	Device	/	25T/Box	CE
		GCDEN(NS)-304a	Card	/		
	WB/S/P	GCDEN(NS)-402a	Device	/	25T/Box	CE
		GCDEN(NS)-402b	Card	/		
		GCDEN(NS)-402c	Card	/		
HBsAb Hepatitis B Rapid Test	S/P	GCHBsb-301a	Strip	30mIU/ml	50T/Box	
		GCHBsb-302a	Device	30mIU/ml	25T/Box	
		GCHBsb-304a	Card	30mIU/ml		
	WB/S/P	GCHBsb-401a	Strip	30mIU/ml	25T/Box	
		GCHBsb-402a	Device	30mIU/ml	50T/Box	
		GCHBsb-404a	Card	30mIU/ml		
HBsAg Hepatitis B Surface Antigen Rapid Test	S/P	GCHBsg-301a	Strip	1ng/ml	100T/Box	CE
		GCHBsg-301b	Dipstick	1ng/ml	100T/Box	CE
		GCHBsg-302a	Device	1ng/ml	25T/Box	CE
		GCHBsg-302b	Device	1ng/ml	25T/Box	CE
		GCHBsg-304a	Card	1ng/ml		
	WB/S/P	GCHBsg-401a	Strip	1ng/ml	50T/Box	CE
		GCHBsg-402a	Device	1ng/ml	25T/Box	CE
		GCHBsg-404a	Card	1ng/ml		
HCV Hepatitis C Virus Rapid Test	S/P	GCHCV-301a	Strip	/	50T/Box	CE
		GCHCV-302a	Device	/	25T/Box	CE
		GCHCV-304a	Card	/		
	WB/S/P	GCHCV-401a	Strip	/	50T/Box	CE
		GCHCV-402a	Device	/	25T/Box	CE
		GCHCV-404a	Card	/		
HIV 1/2 Human Immunodeficiency Virus Rapid Test	S/P	GCHIV-301a	Strip	/	50T/Box	CE
		GCHIV-302a	Device	/	25T/Box	CE
		GCHIV-304a	Card	/		
	WB/S/P	GCHIV-401a	Strip	/	50T/Box	CE
		GCHIV-402a	Device	/	25T/Box	CE
		GCHIV-404a	Card	/		
HIV 1/2/0 Human Immunodeficiency Virus Rapid Test	WB/S/P	GCHIV-T401a	Strip	/	50T/Box	
		GCHIV-T402a	Device	/	25T/Box	
H.pylori Rapid Test	S/P	GCHP-301a	Strip	/	50T/Box	CE
		GCHP-302a	Device	/	25T/Box	CE
		GCHP-304a	Card	/		
	WB/S/P	GCHP-401a	Strip	/	50T/Box	CE
		GCHP-402a	Device	/	25T/Box	CE
		GCHP-404a	Card	/		



Product Name	Specimen	Catalog No.	Format	Sensitivity	Package	Qualification
H.pylori Rapid Test	Feces	GCHP-601a	Strip	/	50T/Box	CE
		GCHP-602a	Device	/	20T/Box	CE
		GCHP-604a	Card	/		
Malaria Pf/Pv Rapid Test	WB	GCMAL(pf/pv)-402a	Device	200 parasites	25T/Box	CE
		GCMAL(pf/pv)-404a	Card	200 parasites		
Malaria Pf/Pan Rapid Test	WB	GCMAL(pf/pan)-402a	Device	200 parasites	25T/Box	CE
		GCMAL(pf/pan)-404a	Card	200 parasites		
Malaria Pf Rapid Test	WB	GCMAL(pf)-401a	Strip	200 parasites	50T/Box	CE
		GCMAL(pf)-402a	Device	200 parasites	25T/Box	CE
		GCMAL(pf)-404a	Card	200 parasites		
Malaria Pan Rapid Test	WB	GCMAL(pan)-401a	Strip	200 parasites	50T/Box	CE
		GCMAL(pan)-402a	Device	200 parasites	25T/Box	CE
		GCMAL(pan)-404a	Card	200 parasites		
Malaria Pf/Pv Ab Rapid Test	S/P	GCMAL(pf/pv Ab)-302a	Device		25T/Box	CE
	WB/S/P	GCMAL(pf/pv Ab)-402a	Device		25T/Box	CE
Gonorrhoea Rapid Test	Swab/Urine	GCGON-501a	Strip	1.0E*7	100T/Box	CE
		GCGON-502a	Device	1.0E*7	25T/Box	CE
Strep A Rapid Test	Throat Swab	GCSTR-501a	Strip	/	50T/Box	CE
		GCSTR-502a	Device	/	20T/Box	CE
		GCSTR-504a	Card	/		
Syphilis Rapid Test	S/P	GCSYP-301a	Strip	/	100T/Box	CE
		GCSYP-301b	Dipstick	/	100T/Box	CE
		GCSYP-302a	Device	/	25T/Box	CE
		GCSYP-302b	Device	/	25T/Box	CE
	WB/S/P	GCSYP-304a	Card	/		
		GCSYP-401a	Strip	/	100T/Box	CE
		GCSYP-401b	Strip	/	100T/Box	CE
		GCSYP-402a	Device	/	25T/Box	CE
		GCSYP-402b	Device	/	25T/Box	CE
		GCSYP-404a	Card	/		
Tuberculosis IgG/IgM Antibody Rapid Test (three lines)	S/P	GCTB-302a	Device	/	25T/Box	CE
		GCTB-304a	Card	/		
	WB/S/P	GCTB-402a	Device	/	25T/Box	CE
		GCTB-404a	Card	/		
Toxo IgG/IgM Rapid Test (three lines)	S/P	GCTOX-302a	Device	/	25T/Box	CE
		GCTOX-304a	Card	/		
	WB/S/P	GCTOX-402a	Device	/	25T/Box	CE
		GCTOX-404a	Card	/		
TOXO IgG Rapid Test	WB/S/P	GCTOX(IgG)-402a	Device	/	25T/Box	CE
		GCTOX(IgG)-404a	Card	/		
TOXO IgM Rapid Test	WB/S/P	GCTOXO(IgM)-402a	Device	/	25T/Box	CE
		GCTOXO(IgM)-404a	Card	/		
Typhoid(Ab) Rapid Test	S/P	GCTYP-302a	Device	/	25T/Box	CE
		GCTYP-304a	Card	/		
S.Typhi Ag Rapid Test	Feces	GCSAL(ST)-602a	Device	/	20T/Box	CE
S.Typhi/paratyphi A Ag Rapid Test	Feces	GCSAL(ST/PTA)-602a	Device	/	20T/Box	CE
CMV IgG Rapid Test	WB/S/P	GCCMV(IgG)-402a	device	/	25T/Box	CE
		GCCMV(IgG)-404a	Card	/		
CMV IgM Rapid Test	S/P	GCCMV(IgM)-302a	device	/	25T/Box	CE
		GCCMV(IgM)-304a	Card	/		
	WB/S/P	GCCMV(IgM)-402a	device	/	25T/Box	CE
		GCCMV(IgM)-404a	Card	/		
HBeAb Rapid Test	S/P	GCHBeb-302a	device	2NCU	25T/Box	
		GCHBeb-302b	device	8NCU	25T/Box	
		GCHBeb-304a	Card	2NCU		
		GCHBeb-304b	Card	8NCU		
HBeAg Rapid Test	S/P	GCHBeg-302a	device	0.5NCU	25T/Box	
		GCHBeg-304a	Card	0.5NCU		



Product Name	Specimen	Catalog No.	Format	Sensitivity	Package	Qualification
HBcAb Rapid Test	S/P	GCHBcb-302a	device	2NCU	25T/Box	
		GCHBcb-302b	device	8NCU	25T/Box	
		GCHBcb-304a	Card	2NCU		
		GCHBcb-304b	Card	8NCU		
HSV1/2 IgG Rapid Test	WB/S/P	GCHSV(IgG)-402a	device	/	25T/Box	CE
		GCHSV(IgG)-404a	Card	/		
HSV1/2 IgM Rapid Test	S/P	GCHSV(IgM)-302a	device	/	25T/Box	CE
		GCHSV(IgM)-304a	Card	/		
	WB/S/P	GCHSV(IgM)-402a	device	/	25T/Box	CE
		GCHSV(IgM)-404a	Card	/		
HSV1 IgG Rapid Test	WB/S/P	GCHSV 1(IgG)-402a	device	/	25T/Box	CE
		GCHSV 1(IgG)-404a	Card	/		
HSV1 IgM Rapid Test	WB/S/P	GCHSV 1(IgM)-402a	device	/	25T/Box	CE
		GCHSV 1(IgM)-404a	Card	/		
HSV2 IgM Rapid Test	WB/S/P	GCHSV 2(IgM)-402a	device	/	25T/Box	CE
		GCHSV 2(IgM)-404a	Card	/		
HSV2 IgG Rapid Test	WB/S/P	GCHSV 2(IgG)-402a	device	/	25T/Box	CE
		GCHSV 2(IgG)-404a	Card	/		
Rota Virus/Adeno Virus Rapid Test	Feces	GCROA/ADE-602a	Device	/	20T/Box	CE
		GCROA/ADE-604a	Card	/		
Rota Virus Rapid Test	Feces	GCROA-602a	device	/	20T/Box	CE
		GCROA-604a	Card	/		
Adeno Virus Rapid Test	Feces	GCADE-602a	device	/	20T/Box	CE
		GCADE-604a	Card	/		
Rubella IgG Rapid Test	WB/S/P	GCRUB(IgG)-402a	device	/	25T/Box	CE
		GCRUB(IgG)-404a	Card	/		
Rubella IgM Rapid Test	S/P	GCRUB(IgM)-302a	Device	/	25T/Box	CE
		GCRUB(IgM)-304a	Card	/		
	WB/S/P	GCRUB(IgM)-402a	device	/	25T/Box	CE
		GCRUB(IgM)-404a	Card	/		
HAV Rapid Test	S/P	GCHAV-302a	Device	/	25T/Box	CE
HAV Ag Rapid Test	Feces	GCHAV-602a	Device	/	20T/Box	CE
HEV IgM Rapid Test	S/P	GCHEV-302a	Device	/	25T/Box	CE
Flu A Ag Rapid Test	Nasal /	GCFLU(A)-501a	Strip	1.5 x 10 ⁴ TCID ₅₀		CE
	Throat Swabs	GCFLU(A)-502a	Device	1.5 x 10 ⁴ TCID ₅₀	15T/Box	CE
Flu A/B Ag Rapid Test	Nasal / Throat Swabs	GCFLU(A/B)-501a	Strip	1.5x 10 ⁴ TCID ₅₀ / 1.5 x 10 ⁵ TCID ₅₀		CE
		GCFLU(A/B)-502a	Device	1.5x 10 ⁴ TCID ₅₀ / 1.5 x 10 ⁵ TCID ₅₀	15T/Box	CE
V.cholerae O1 Ag Rapid Test	Feces	GEVCH(O1)-602a	Device	/	20T/Box	CE
V.cholerae O1/O139 Ag Rapid Test	Feces	GEVCH(O1/O9)-602a	Device	/	20T/Box	CE
Chikungunya IgM Rapid Test	S/P	GCCHK(IgM)-302a	Device	/	25T/Box	CE
	WB/S/P	GCCHK(IgM)-402a	Device	/	25T/Box	CE
Chagas Ab Rapid Test	S/P	GCCHA-302a	Device	/	25T/Box	CE
	WB/S/P	GCCHA-402a	Device	/	25T/Box	CE

Multi-ID Device

Product Name	Specimen	Catalog No.	Combination	Package
HBV Combo Rapid Test	S/P	GCHBV-355	GCHBsg, GCHBsb, GCHBeg, GCHBeb, GCHBcb	20T/Box
TORCH IgM Combo Rapid Test	S/P	GCTOM-345a	Toxo IgM, Rubella IgM, CMV IgM, HSV IgM	20T/Box
TORCH IgG Combo Rapid Test	S/P	GCTOG-345a	Toxo IgG, Rubella IgG, CMV IgG, HSV IgG	20T/Box
Dengue NS1-IgG/IgM Combo Rapid Test	WB/S/P	GC DEN-425a	Dengue NS1, Dengue IgG/IgM	20T/Box

Drugs of Abuse



Product Name	Specimen	Catalog No.	Format	Sensitivity	Package	Qualification
AMP (Amphetamine) Rapid Test	Urine	GBAMP-101a	Strip	1000ng/ml	100T/Box	CE
		GBAMP-101b	Strip	2000ng/ml	100T/Box	CE
		GBAMP-U101a	Strip	500ng/ml	100T/Box	CE
		GBAMP-A101a	Strip	300ng/ml	100T/Box	CE
		GBAMP-102a	Device	1000ng/ml	25T/Box	CE
		GBAMP-102b	Device	2000ng/ml	25T/Box	CE
		GBAMP-A102a	Device	300ng/ml	25T/Box	CE
		GBAMP-U102a	Device	500ng/ml	25T/Box	CE
		GBAMP-104a	Card	1000ng/ml	25T/Box	
		GBAMP-106a	Cup	1000ng/ml	25T/Box	CE
		GBAMP-104b	Card	2000ng/ml		
		GBAMP-106b	Cup	2000ng/ml	25T/Box	CE
		GBAMP-U104a	Card	500ng/ml		
		GBAMP-U106a	Cup	500ng/ml	25T/Box	CE
		GBAMP-U104a	Card	500ng/ml		
		GBAMP-A104a	Card	300ng/ml		
		GBAMP-A106a	Cup	300ng/ml	25T/Box	CE
		GBAMP-105a	Panel	1000ng/ml	25T/Box	CE
		GBAMP-105b	Panel	2000ng/ml	25T/Box	CE
		GBAMP-U105a	Panel	500ng/ml	25T/Box	CE
GBAMP-A105a	Panel	300ng/ml	25T/Box	CE		
BAR (Barbiturates) Rapid Test	Urine	GBBAR-101a	Strip	300ng/ml	100T/Box	CE
		GBBAR-101b	Strip	600ng/ml	100T/Box	CE
		GBBAR-A101a	Strip	200ng/ml	100T/Box	CE
		GBBAR-102a	Device	300ng/ml	25T/Box	CE
		GBBAR-102b	Device	600ng/ml	25T/Box	CE
		GBBAR-A102a	Device	200ng/ml	25T/Box	CE
		GBBAR-104a	Card	300ng/ml		
		GBBAR-106a	Cup	300ng/ml	25T/Box	CE
		GBBAR-104b	Card	600ng/ml		
		GBBAR-106b	Cup	600ng/ml	25T/Box	CE
		GBBAR-A104a	Card	200ng/ml		
		GBBAR-A106a	Cup	200ng/ml	25T/Box	CE
		GBBAR-105a	Panel	300ng/ml	25T/Box	CE
		GBBAR-105b	Panel	600ng/ml	25T/Box	CE
GBBAR-A105a	Panel	200ng/ml	25T/Box	CE		
BUP (Buprenorphine) Rapid Test	Urine	GBBUP-101a	Strip	10ng/ml	100T/Box	CE
		GBBUP-102a	Device	10ng/ml	25T/Box	CE
		GBBUP-104a	Card	10ng/ml		
		GBBUP-106a	Cup	10ng/ml	25T/Box	CE
		GBBUP-105a	Panel	10ng/ml	25T/Box	CE
BZO (Benzodiazepines) Rapid Test	Urine	GBBZO-101a	Strip	300ng/ml	100T/Box	CE
		GBBZO-101b	Strip	600ng/ml	100T/Box	CE
		GBBZO-A101a	Strip	200ng/ml	100T/Box	CE
		GBBZO-102a	Device	300ng/ml	25T/Box	CE
		GBBZO-102b	Device	600ng/ml	25T/Box	CE
		GBBZO-A102a	Device	200ng/ml	25T/Box	CE
		GBBZO-104b	Card	600ng/ml		
		GBBZO-104a	Card	300ng/ml		



Product Name	Specimen	Catalog No.	Format	Sensitivity	Package	Qualification
BZO (Benzodiazepines) Rapid Test	Urine	GBBZO-106a	Cup	300ng/ml	25T/Box	CE
		GBBZO-106b	Cup	600ng/ml	25T/Box	CE
		GBBZO-A104a	Card	200ng/ml		
		GBBZO-A106a	Cup	200ng/ml	25T/Box	CE
		GBBZO-105a	Panel	300ng/ml	25T/Box	CE
		GBBZO-105b	Panel	600ng/ml	25T/Box	CE
		GBBZO-A105a	Panel	200ng/ml	25T/Box	CE
COC (Cocaine) Rapid Test	Urine	GBCOC-101a	Strip	300ng/ml	100T/Box	CE
		GBCOC-101b	Strip	600ng/ml	100T/Box	CE
		GBCOC-U101a	Strip	150ng/ml	100T/Box	CE
		GBCOC-102a	Device	300ng/ml	25T/Box	CE
		GBCOC-102b	Device	600ng/ml	25T/Box	CE
		GBCOC-104a	Device	300ng/ml	25T/Box	CE
		GBCOC-U102a	Device	150ng/ml	25T/Box	CE
		GBCOC-106a	Cup	300ng/ml	25T/Box	CE
		GBCOC-104b	Card	600ng/ml		
		GBCOC-106b	Cup	600ng/ml	25T/Box	CE
		GBCOC-U104a	Card	150ng/ml		
		GBCOC-U106a	Cup	150ng/ml	25T/Box	CE
		GBCOC-105a	Panel	300ng/ml	25T/Box	CE
		GBCOC-105b	Panel	600ng/ml	25T/Box	CE
GBCOC-U105a	Panel	150ng/ml	25T/Box	CE		
COT (Cotinine) Rapid Test	Urine	GBCOT-102a	Device	200ng/ml	25T/Box	CE
		GBCOT-104a	Card	200ng/ml		
		GBCOT-106a	Cup	200ng/ml	25T/Box	CE
		GBCOT-105a	Panel	200ng/ml	25T/Box	CE
K2 (Synthetic Marijuana) Rapid Test	Urine	GBK2-101a	Strip	50ng/ml	100T/Box	CE
		GBK2-102a	Device	50ng/ml	25T/Box	CE
		GBK2-101b	Strip	25ng/ml	100T/Box	CE
		GBK2-102b	Strip	25ng/ml	100T/Box	CE
		GBK2-104a	Card	50ng/ml		
		GBK2-106a	Cup	50ng/ml	25T/Box	CE
		GBK2-104b	Card	25ng/ml		
		GBK2-106b	Cup	25ng/ml	25T/Box	CE
		GBK2-105a	Panel	50ng/ml	25T/Box	CE
GBK2-105b	Panel	25ng/ml	25T/Box	CE		
KET (Ketamine) Rapid Test	Urine	GBKET-101a	Strip	1000ng/ml	100T/Box	CE
		GBKET-101b	Strip	2000ng/ml	100T/Box	CE
		GBKET-101c	Strip	3000ng/ml	100T/Box	CE
		GBKET-101d	Strip	100ng/ml	100T/Box	CE
		GBKET-102a	Device	1000ng/ml	25T/Box	CE
		GBKET-102b	Device	2000ng/ml	25T/Box	CE
		GBKET-102c	Device	3000ng/ml	25T/Box	CE
		GBKET-102d	Device	100ng/ml	25T/Box	CE
		GBKET-104a	Card	1000ng/ml		
		GBKET-106a	Cup	1000ng/ml	25T/Box	CE
		GBKET-104a	Card	1000ng/ml		
		GBKET-104b	Card	2000ng/ml		
		GBKET-106b	Cup	2000ng/ml	25T/Box	CE
		GBKET-104c	Card	3000ng/ml		
		GBKET-106c	Cup	3000ng/ml	25T/Box	CE
		GBKET-104d	Card	100ng/ml		
		GBKET-106d	Cup	100ng/ml	25T/Box	CE
GBKET-105a	Panel	1000ng/ml	25T/Box	CE		
GBKET-105b	Panel	2000ng/ml	25T/Box	CE		
GBKET-105c	Panel	3000ng/ml	25T/Box	CE		
GBKET-105d	Panel	100ng/ml	25T/Box	CE		
MDMA (Ecstasy) Rapid Test	Urine	GBMDM-101a	Strip	500ng/ml	100T/Box	CE
		GBMDM-101b	Strip	2000ng/ml	100T/Box	CE



Drugs of Abuse

Product Name	Specimen	Catalog No.	Format	Sensitivity	Package	Qualification
MDMA (Ecstasy) Test	Urine	GBMDM-101c	Strip	1000ng/ml	100T/Box	CE
		GBMDM-102a	Device	500ng/ml	25T/Box	CE
		GBMDM-102b	Device	2000ng/ml	25T/Box	CE
		GBMDM-102c	Device	1000ng/ml	25T/Box	CE
		GBMDM-104a	Card	500ng/ml		
		GBMDM-106a	Cup	500ng/ml	25T/Box	CE
		GBMDM-104b	Card	2000ng/ml		
		GBMDM-106b	Cup	2000ng/ml	25T/Box	CE
		GBMDM-104c	Card	1000ng/ml		
		GBMDM-106c	Cup	1000ng/ml	25T/Box	CE
		GBMDM-105a	Strip	500ng/ml	100T/Box	CE
		GBMDM-105b	Strip	2000ng/ml	100T/Box	CE
		GBMDM-105c	Strip	1000ng/ml	100T/Box	CE
		MET (Methamphetamine) Rapid Test	Urine	GBMET-101a	Strip	1000ng/ml
GBMET-101b	Strip			2000ng/ml	100T/Box	CE
GBMET-U101a	Strip			500ng/ml	100T/Box	CE
GBMET-A101a	Strip			300ng/ml	100T/Box	CE
GBMET-102a	Device			1000ng/ml	25T/Box	CE FDA
GBMET-102b	Device			2000ng/ml	25T/Box	CE
GBMET-U102a	Device			500ng/ml	25T/Box	CE
GBMET-A102a	Device			300ng/ml	25T/Box	CE
GBMET-104a	Card			1000ng/ml		
GBMET-106a	Cup			1000ng/ml	25T/Box	CE FDA
GBMET-104b	Card			2000ng/ml		
GBMET-106b	Cup			2000ng/ml	25T/Box	CE
GBMET-U104a	Card			500ng/ml		
GBMET-U106a	Cup			500ng/ml	25T/Box	CE
GBMET-A104a	Card			300ng/ml		
GBMET-A106a	Cup			300ng/ml	25T/Box	CE
GBMET-105a	Panel			1000ng/ml	25T/Box	CE FDA
GBMET-105b	Panel			2000ng/ml	25T/Box	CE
GBMET-U105a	Panel	500ng/ml	25T/Box	CE		
GBMET-A105a	Panel	300ng/ml	25T/Box	CE		
MOP (Morphine) Rapid Test	Urine	GBMOP-101a	Strip	300ng/ml	100T/Box	CE
		GBMOP-101b	Strip	600ng/ml	100T/Box	CE
		GBMOP-101c	Device	2000ng/ml	25T/Box	CE
		GBMOP-102a	Device	300ng/ml	25T/Box	CE
		GBMOP-102b	Device	600ng/ml	25T/Box	CE
		GBMOP-102c	Device	2000ng/ml	25T/Box	CE
		GBMOP-104a	Card	300ng/ml		
		GBMOP-106a	Cup	300ng/ml	25T/Box	CE
		GBMOP-104b	Card	600ng/ml		
		GBMOP-106b	Cup	600ng/ml	25T/Box	CE
		GBMOP-104c	Card	2000ng/ml		
		GBMOP-106c	Cup	2000ng/ml	25T/Box	CE
		GBMOP-105a	Panel	300ng/ml	25T/Box	CE
		GBMOP-105b	Panel	600ng/ml	25T/Box	CE
GBMOP-105c	Panel	2000ng/ml	25T/Box	CE		
MLQ (Methaqualone) Rapid Test	Urine	GBMQL-101a	Strip	300ng/ml	100T/Box	CE
		GBMQL-102a	Device	300ng/ml	25T/Box	CE
		GBMQL-104a	Card	300ng/ml		
		GBMQL-106a	Cup	300ng/ml	25T/Box	CE
		GBMQL-105a	Panel	300ng/ml	25T/Box	CE
MTD (Methadone) Rapid Test	Urine	GBMTD-101a	Strip	300ng/ml	100T/Box	CE
		GBMTD-A101a	Strip	200ng/ml	100T/Box	CE
		GBMTD-101b	Strip	600ng/ml	100T/Box	CE
		GBMTD-102a	Device	300ng/ml	25T/Box	CE
		GBMTD-A102a	Device	200ng/ml	25T/Box	CE
		GBMTD-102b	Device	600ng/ml	25T/Box	CE



Product Name	Specimen	Catalog No.	Format	Sensitivity	Package	Qualification
MTD (Methadone) Rapid Test	Urine	GBMTD-104a	Card	300ng/ml		
		GBMTD-106a	Cup	300ng/ml	25T/Box	CE
		GBMTD-A104a	Card	200ng/ml		
		GBMTD-A106a	Cup	200ng/ml	25T/Box	CE
		GBMTD-104b	Card	600ng/ml		
		GBMTD-106b	Cup	600ng/ml	25T/Box	CE
		GBMTD-105a	Panel	300ng/ml	25T/Box	CE
		GBMTD-A105a	Panel	200ng/ml	25T/Box	CE
		GBMTD-105b	Panel	600ng/ml	25T/Box	CE
OPI (Opiate) Rapid Test	Urine	GBOPI-101a	Strip	2000ng/ml	100T/Box	CE
		GBOPI-102a	Device	2000ng/ml	25T/Box	CE
		GBOPI-104a	Card	2000ng/ml		
		GBOPI-106a	Cup	2000ng/ml	25T/Box	CE
		GBOPI-105a	Panel	2000ng/ml	25T/Box	CE
OXY (Oxycodone) Rapid Test	Urine	GBOXY-101a	Strip	100ng/ml	100T/Box	CE
		GBOXY-102a	Device	100ng/ml	25T/Box	CE
		GBOXY-104a	Card	100ng/ml		
		GBOXY-106a	Cup	100ng/ml	25T/Box	CE
		GBOXY-105a	Panel	500ng/ml	25T/Box	CE
PCP (Phencyclidine) Rapid Test	Urine	GBPCP-101a	Strip	25ng/ml	100T/Box	CE
		GBPCP-104a	Device	25ng/ml	25T/Box	CE
		GBPCP-102a	Device	25ng/ml	25T/Box	CE
		GBPCP-101b	Strip	50ng/ml	100T/Box	CE
		GBPCP-102b	Device	50ng/ml	25T/Box	CE
		GBPCP-104a	Card	25ng/ml		
		GBPCP-106a	Cup	25ng/ml	25T/Box	CE
		GBPCP-104b	Card	50ng/ml		
		GBPCP-106b	Cup	50ng/ml	25T/Box	CE
		GBPCP-105a	Panel	25ng/ml	25T/Box	CE
GBPCP-105b	Panel	50ng/ml	25T/Box	CE		
PPX (Propoxyphene) Rapid Test	Urine	GBPPX-101a	Strip	300ng/ml	100T/Box	CE
		GBPPX-101b	Strip	600ng/ml	100T/Box	CE
		GBPPX-102a	Device	300ng/ml	25T/Box	CE
		GBPPX-102b	Device	600ng/ml	25T/Box	CE
		GBPPX-104a	Card	300ng/ml		
		GBPPX-106a	Cup	300ng/ml	25T/Box	CE
		GBPPX-104b	Card	600ng/ml		
		GBPPX-106b	Cup	600ng/ml	25T/Box	CE
		GBPPX-105a	Panel	300ng/ml	25T/Box	CE
		GBPPX-105b	Panel	600ng/ml	25T/Box	CE
TCA (Tricyclic Antidepressants) Rapid Test	Urine	GBTCA-101a	Strip	1000ng/ml	100T/Box	CE
		GBTCA-102a	Device	1000ng/ml	25T/Box	CE
		GBTCA-104a	Card	1000ng/ml		
		GBTCA-106a	Cup	1000ng/ml	25T/Box	CE
		GBTCA-105a	Panel	1000ng/ml	25T/Box	CE
THC (Marijuana) Rapid Test	Urine	GBTHC-101a	Strip	50ng/ml	100T/Box	CE FDA
		GBTHC-101b	Strip	200ng/ml	100T/Box	CE
		GBTHC-101c	Strip	20ng/ml	100T/Box	CE
		GBTHC-A101a	Strip	15ng/ml	100T/Box	CE
		GBTHC-B101a	Strip	150ng/ml	100T/Box	CE
		GBTHC-102a	Device	50ng/ml	25T/Box	CE FDA
		GBTHC-102b	Device	200ng/ml	25T/Box	CE
		GBTHC-102c	Device	20ng/ml	25T/Box	CE
		GBTHC-A102a	Device	15ng/ml	25T/Box	CE
		GBTHC-B102a	Device	150ng/ml	25T/Box	CE
		GBTHC-104a	Card	50ng/ml		
		GBTHC-106a	Cup	50ng/ml	25T/Box	CE FDA
		GBTHC-104b	Card	200ng/ml		
		GBTHC-106b	Cup	200ng/ml	25T/Box	CE



Drugs of Abuse

Product Name	Specimen	Catalog No.	Format	Sensitivity	Package	Qualification
THC (Marijuana) Rapid Test	Urine	GBTHC-104c	Card	20ng/ml		
		GBTHC-106c	Cup	20ng/ml	25T/Box	CE
		GBTHC-A104a	Card	15ng/ml		
		GBTHC-A106a	Cup	15ng/ml	25T/Box	CE
		GBTHC-A104a	Card	15ng/ml		
		GBTHC-B104a	Card	150ng/ml		
		GBTHC-B106a	Cup	150ng/ml	25T/Box	CE
		GBTHC-105a	Panel	50ng/ml	25T/Box	CE FDA
		GBTHC-105b	Panel	200ng/ml	25T/Box	CE
		GBTHC-105c	Panel	20ng/ml	25T/Box	CE
		GBTHC-A105a	Panel	15ng/ml	25T/Box	CE
		GBTHC-B105a	Panel	150ng/ml	25T/Box	CE
		TRA (Tramadol) Rapid Test	Urine	GBTRA-101a	Strip	200ng/ml
GBTRA-102a	Device			200ng/ml	25T/Box	CE
GBTRA-104a	Card			200ng/ml		
GBTRA-106a	Cup			200ng/ml	25T/Box	CE
GBTRA-105a	Panel			200ng/ml	25T/Box	CE
FEN (Fentanyl) Rapid Test	Urine	GBFEN-101a	Strip	300ng/ml	100T/Box	CE
		GBFEN-102a	Device	300ng/ml	25T/Box	CE
		GBFEN-104a	Card	300ng/ml		
		GBFEN-106a	Cup	300ng/ml	25T/Box	CE
		GBFEN-105a	Panel	300ng/ml	25T/Box	CE
EDDP (Methadone Metabolite) Rapid Test	Urine	GBEDD-101a	Strip	300ng/ml	100T/Box	CE
		GBEDD-101b	Strip	100ng/ml	100T/Box	CE
		GBEDD-102a	Device	300ng/ml	25T/Box	CE
		GBEDD-102b	Device	100ng/ml	25T/Box	CE
		GBEDD-104a	Card	300ng/ml		
		GBEDD-104b	Card	100ng/ml		
		GBEDD-106a	Cup	300ng/ml	25T/Box	CE
		GBEDD-106b	Cup	100ng/ml	25T/Box	CE
		GBEDD-105a	Panel	300ng/ml	25T/Box	CE
GBEDD-105b	Panel	100ng/ml	25T/Box	CE		
Saliva Alcohol Test Strip	Urine	DAL-701a	Strip	0.02%	100T/Box	

Drugs of Abuse Panel

Combination	Specimen	Inside Catalog No.	Format	Package	Qualification
MOP/MDMA/BUP	Urine	GBDOA-134-005	Panel	25T/Box	CE
MOP/BZO/THC	Urine	GBDOA-134-008	Panel	25T/Box	CE
MOP/EDDP/BUP	Urine	GBDOA-134-010	Panel	25T/Box	CE
AMP/MET/MDMA/COC/THC/BAR/BZO/MOP/OPI/MTD	Urine	GBDOA-1104-003	Panel	20T/Box	CE
AMP/MET/MDMA/COC/THC/BAR/BZO/MOP/BUP/MTD	Urine	GBDOA-1104-009	Panel	20T/Box	CE
AMP/MET/MDMA/COC/THC/TRM/BZO/MOP/BUP/MTD	Urine	GBDOA-1104-013	Panel	20T/Box	CE
AMP1000ng/ml,THC200ng/ml,MOP 600ng/ml	Urine	GBDOA-134-001	Panel	25T/Box	CE
AMP1000ng/ml,THC200ng/ml,MOP600ng/ml, MET2000ng/ml,COC600ng/ml,BZO,600ng/ml	Urine	GBDOA164-001	Panel	20T/Box	CE
AMP1000ng/ml,BAR,BZO,COC,THC,MTD,MET,MDMA, MOP,TCA	Urine	GBD1104-001	Panel	20T/Box	CE
AMP1000ng/ml,BAR,BZO,COC,THC,MTD,MET,MDMA, MOP,TCA,OPI	Urine	GBDOA1124-001	Panel	20T/Box	CE
MOP300,MET1000ng/ml	Urine	GBDOA124-001	Panel	25T/Box	CE



Rapid Test



Pregnancy & Fertility

Product Name	Specimen	Catalog No.	Format	Sensitivity	Package	Qualification
hCG Pregnancy Rapid Test	Urine	GAHCG-101a	Strip (Dipstick)	25mIU/ml	100T/Box	CE 0197 FDA
		GAHCG-102a	Device	25mIU/ml	25T/Box	CE 0197 FDA
		GAHCG-103a	Midstream	25mIU/ml	20T/Box	CE 0197 FDA
		GAHCG-104a	Card	25mIU/ml		
		GAHCG-101d	Strip (Dipstick)	20mIU/ml	100T/Box	CE
		GAHCG-102d	Device	20mIU/ml	25T/Box	CE
		GAHCG-103d	Midstream	20mIU/ml	20T/Box	CE
		GAHCG-101e	Strip (Dipstick)	25mIU/ml, Detect Pregnancy five days before next Period	100T/Box	CE 0197
		GAHCG-102e	Device	25mIU/ml, Detect Pregnancy five days before next Period	25T/Box	CE 0197
		GAHCG-103e	Midstream	25mIU/ml, Detect Pregnancy five days before next Period	20T/Box	CE 0197
		GAHCG-104d	Card	25mIU/ml		
		GAHCG-201a	Strip	25mIU/ml	100T/Box	CE
		GAHCG-202a	Device	25mIU/ml	25T/Box	CE
		GAHCG-201b	Strip	10mIU/ml	100T/Box	CE
GAHCG-202b	Device	10mIU/ml	25T/Box	CE		
GAHCG-204b	Card	10mIU/ml				
GAHCG-204a	Card	25mIU/ml				
hCG Ultra Pregnancy Rapid Test	Urine	GAHCG-101b	Strip(Dipstick)	10mIU/ml	100T/Box	CE
		GAHCG-102b	Device	10mIU/ml	25T/Box	CE
		GAHCG-103b	Midstream	10mIU/ml	20T/Box	CE
		GAHCG-104b	Card	10mIU/ml		
LH Ovulation Rapid Test	Urine	GALH-101a	Strip(Dipstick)	25mIU/ml	100T/Box	CE 0197 FDA
		GALH-102a	Device	25mIU/ml	25T/Box	CE 0197 FDA
		GALH-103a	Midstream	25mIU/ml	20T/Box	CE 0197 FDA
		GALH-104a	Card	25mIU/ml		
		GALH-101b	Strip(Dipstick)	40mIU/ml	100T/Box	CE FDA
		GALH-102b	Device	40mIU/ml	25T/Box	CE FDA
		GALH-103b	Midstream	40mIU/ml	20T/Box	FDA
		GALH-101c	Strip	20mIU/ml		CE FDA
		GALH-103c	Midstream	20mIU/ml		FDA
FSH Menopause Rapid Test	Urine	GAFSH-101a	Strip	25mIU/ml	100T/Box	CE
		GAFSH-102a	Device	25mIU/ml	25T/Box	CE
IGFBP-1 PROM Rapid Test	Cervical Secretion	GAIGF1-101a	Strip	25ng/ml	100T/Box	
		GAIGF1-102a	Device	25ng/ml	25T/Box	

Rapid Test



Tumor Marker

Product Name	Specimen	Catalog No.	Format	Sensitivity	Package	Qualification
FOB Fecal Occult Blood Rapid Test	Feces	GEFOB-601a	Strip	50ng/ml	50T/Box	CE FDA
		GEFOB-601b	Strip	100ng/ml	50T/Box	CE
		GEFOB-601c	Strip	50ng/ml	50T/Box	CE
		GEFOB-604a	Card	50ng/ml		
		GEFOB-602a	Device	50ng/ml	25T/Box	CE FDA
		GEFOB-602b	Device	100ng/ml	25T/Box	CE
PSA Prostate Specific Antigen Rapid Test	S/P	GEPSA-302a	Device	4ng/ml	25T/Box	CE 0197
		GEPSA-301a	Strip	4ng/ml	100T/Box	CE 0197
		GEPSA-304a	Card	4ng/ml		
	WB/S/P	GEPSA-402a	Device	4ng/ml	25T/Box	CE
		GEPSA-401a	Strip	4ng/ml	100T/Box	CE
		GEPSA-404a	Card	4ng/ml		
PSA Prostate Specific Antigen Semi-Quantitative Rapid Test	WB/S/P	GEPSA-402b	Device	4ng/ml,10ng/ml	50T/Box	
		GEPSA-404b	Card	4ng/ml,10ng/ml		
Alpha Fetal Protein Rapid Test	S/P	GEAFP-302a	Device	20ng/ml	25T/Box	CE
		GEAFP-304a	Card	20ng/ml		
	WB/S/P	GEAFP-402a	Device	20ng/ml	25T/Box	CE
		GEAFP-404a	Card	20ng/ml		
CEA Carcinoembryonic Antigen Rapid Test	S/P	GECEA-302a	Device	5ng/ml	25T/Box	CE
		GECEA-304a	Card	5ng/ml		
	WB/S/P	GECEA-402a	Device	5ng/ml	25T/Box	CE
		GECEA-404a	Card	5ng/ml		
Transferrin (TF) Rapid Test	Feces	GETF-601a	Strip	10ng/ml	50T/Box	CE
		GETF-602a	Device	10ng/ml	20T/Box	CE
FOB / Transferrin (TF) Combo Rapid Test	Feces	GEFOB/TF-602a	Device	50ng/ml/10ng/ml	20T/Box	CE

New and Precise
High Sensitivity and Specificity
Stable Quality





Rapid Test

Cardiac Marker



Product Name	Specimen	Catalog No.	Format	Sensitivity	Package	Qualification
Troponin I Rapid Test	S/P	GDTRO-301a	Strip	0.5ng/ml	50T/Box	CE
		GDTRO-302a	Device	0.5ng/ml	25T/Box	CE
		GDTRO-304a	Card	0.5ng/ml		
	WB/S/P	GDTRO-401a	Strip	0.5ng/ml	50T/Box	CE
		GDTRO-402a	Device	0.5ng/ml	25T/Box	CE
		GDTRO-404a	Card	0.5ng/ml		
Myoglobin Rapid Test	WB/S/P	GDMYO-402a	Device	50ng/ml	25T/Box	CE
CK-MB Rapid Test	S/P	GDCKM-302a	Device	5ng/ml	25T/Box	CE
		GDCKM-304a	Device	5ng/ml	25T/Box	CE
	WB/S/P	GDCKM-402a	Device	5 ng/ml	25T/Box	CE
		GDCKM-404a	Device	5 ng/ml	25T/Box	CE
CRP Semi-Quantitative C-Reactive Protein Rapid Test	WB/S/P	GDCRP-402a	Device	10mg/ml	25T/Box	CE
		GDCRP-404a	Card	10mg/ml		
Cardiac Myoglobin/CK- MB/ cTnI Combo Rapid Test	WB/S/P	GDCAR-435a	Device	50/5/0.5ng/ml	25T/Box	CE

Rapid Test Brand





Urine Analyzer & Urinalysis Strip

Urine Analyzer

1. Testing Principle: super-high luminosity cold light source reflection testing principle
2. Testing Items: Leucocyte (LEU) Nitrite (NIT) Urobilinogen (URO) Protein (PRO) pH Blood (BLD) Specific Gravity (SG) Ascorbate (VC) Ketone (KET) Bilirubin (BIL) Glucose (GLU)
3. Matching Test Strip: URS-10T, URS-11
4. Light Wavelength: 420nm, 525nm, 560nm, 610nm, 660nm, 950nm
5. Testing Speed: 120 times/h for single testing, 500 times/h for continuous testing.
6. Data Memory: 1000 test results
7. Display Screen: W-200B 240mm*128mm, 5.7" LED screen
8. Printer: Thermal printer inside, with outer stylus printer interface
9. External Output: Connect with the computer through Standard RS232C port
10. Can connect with urine sediments analyzer
11. Can automatically modify the influence to the test results caused by test paper's urine-PH, specific gravity and color
12. Baud Rate of Data Communication Interface: 9600bps
13. Language: English
14. Environment Condition 0~40°C: RH<85%
15. Power Supply AC 220V(±15%), 50~60Hz
16. Fuse wire specification: 250V 2A
17. Power: 60W
18. Net: 3Kg




Urinalysis Strip

Brand




Product Name	Specimen	Format	Sensitivity	Reading Time	Qualification
Glucose	Urine	Strip	2.8~5.5mmol/L	30s	CE FDA
Protein	Urine	Strip	0.15~0.3g/L	60s	CE FDA
Ketone	Urine	Strip	0.5~1.0mmol/L	40s	CE FDA
Blood	Urine	Strip	5-15Ery/μL	60s	CE FDA
Leukocytes	Urine	Strip	5-15Leuko/μL	120s	CE FDA
Urobilinogen	Urine	Strip	3.3-16μmol/L	60s	CE FDA
Bilirubin	Urine	Strip	8.6- 17μmol/L	60s	CE FDA
pH	Urine	Strip	0.5	60s	CE FDA
Specific Gravity	Urine	Strip	0.005	60s	CE FDA
Ascorbate	Urine	Strip	0.5-0.6mmol/L	60s	CE FDA
Nitrite	Urine	Strip	13~22μmol/L	60s	CE FDA
Micro Albumin	Urine	Strip	0.08~0.15mg/dL	60S	CE

Veterinary Test



Product Name	Test Sample	Pack
CCPV Ag Rapid Test	Feces	25test/box
CDV Ag Rapid Test	Serum/Plasma	25test/box
CCV Ag Rapid Test	Feces	25test/box
CCV/CPV Ag Rapid Test	Feces	25test/box
CIV Ag Rapid Test	Nasal Swab	25test/box
CAV Ag Rapid Test	Body secretions	25test/box
FPV Rapid Test	Feces	25test/box
FeLV Rapid Test	WB/Serum/Plasma	25test/box
FIV Rapid Test	WB/Serum/Plasma	25test/box
TgAb Rapid Test	WB/Serum/Plasma	25test/box
Rota Rapid Test	Feces	25test/box
FCoV Rapid Test	Feces	25test/box
ICH Rapid Test	WB/Serum/Plasma	25test/box
CPV Ab Rapid Test	WB/Serum/Plasma	25test/box
CDV Ab Rapid Test	WB/Serum/Plasma	25test/box
ICH Ab Rapid Test	WB/Serum/Plasma	25test/box
cTnl Rapid Test	WB/Serum/Plasma	25test/box
D-dimer Rapid Test	WB/Serum/Plasma	25test/box



Healgen Scientific Limited Liability Company

Add: 3818 Fuqua street, Houston, USA. TX 77047

Tel: +1 713-733-8088

Toll free: 866 982 3818

Fax: +1 713-733-8848

E-mail: sales@healgen.com

ZheJiang Orient Gene Biotech Co., Ltd

Add: 679#, East Yangguang Avenue, Dipu Town, Anji

Tel: +86-572-5303755

Fax: +86-572-5226222

Web: <http://www.orientgene.com>

P.C.: 313300

E-mail: sales@orientgene.com

